

IN-LINE POWER SUPPLIES



AC150A Power Supply



AC137A Power Supply

All AC series in-line power supplies are CE marked for the EC Low Voltage Directive and EMC Directive, and all have UL and TUV approval. The AC150A has standard IEC power input plugs and operate over mains power ratings of 100-240 VAC, 50-60 Hz. AC150A includes a USA, EURO, or Chinese power cord. (ACCORD US/EURO/CN).

- | | | |
|---------------|--------------------------|--|
| AC101A | ±12 volt, +5 volt, 1 amp | Connects the LDF100C to the AC mains wall outlet. One supply is included with each LDF100C module. |
| AC137A | +6 volt, 1.5 amp | Powers the heating element for any of the TSD137 series pneumotachs. The AC137A has fixed built-in USA-style power prongs and does not need a power adapter cord for USA operation. Interchangeable prongs are included for international use. |
| AC150A | +12 volt, 3.4 amp | Connects the MP160/150 System or GASSYSTEM2 to the AC mains wall outlet. One supply is included with each MP160/150 Starter system or GASSYSTEM2. |
| AC300A | +12 volt, 1.25 amp | Connects the MP to mains wall outlet. One supply is included with each MP36/35 system. |

See also: IPS100C Isolated Power Supply

NOTE: The older style AC137A with the ACCORD US/EURO power cord was discontinued in October of 2017.