

## Application Note 236 BIOPAC Software Release History & OS Compatibility

If the information you require isn't listed, please contact BIOPAC at [info@biopac.com](mailto:info@biopac.com) with a specific release question.

### ACQKNOWLEDGE SOFTWARE

#### Windows

Version	Build Date	Supported OS	Hardware
5.0.8.1	10.20.23	10, 11	MP160/150/36R, B-Alert, BHBT, Eye Tracking, Mobita/BioNomadix/BN-Logger/Smart Center/Stellar (real-time and memory implants)
5.0.8	09.13.2022	10, 11	MP160/150/36R, B-Alert, BHBT, Eye Tracking, Mobita/BioNomadix/BN-Logger/Smart Center/Stellar
5.0.7	12.28.2021	10	MP160/150/36R, B-Alert, BHBT, Eye Tracking, Mobita/BioNomadix/BN-Logger/Smart Center/Stellar
5.0.6	06.02.2021	10, 8, 7 (64-bit only)	MP160/150/36R, B-Alert, BHBT, Eye Tracking, Mobita/BioNomadix/BN-Logger/Smart Center/Stellar
5.0.5	08.19.2019	10, 8, 7 (64-bit only)	MP160/150/36R, B-Alert, BHBT, Mobita/BioNomadix/BN-Logger/Smart Center/Stellar
5.0.4	01.11.2019	10, 8, 7 (64-bit only)	MP160/150/36R, B-Alert, BHBT, Mobita/BioNomadix/BN-Logger/Smart Center/Stellar
5.0.3	09.07.2018	10, 8, 7 (64-bit only)	MP160/150/36R, B-Alert, BHBT, Mobita/BioNomadix/BN-Logger/Smart Center/Stellar
5.0.2.1	12.13.2017	10, 8, 7 (64-bit only)	MP160/150/36R, BHBT, Mobita/BioNomadix/BN-Logger/Smart Center/Stellar
5.0.1	11.14.2016	10, 8.x, 7 (64 bit only)	MP160/150/36R, BHBT, Mobita, BioNomadix/BN-Logger/Stellar
5.0	07.28.2016	10, 8.x, 7 (64 bit only)	MP160/150/36R, BHBT, Mobita, BioNomadix/BN-Logger
4.4.2	04.09.2016	10, 8.x & 7	MP150/36R, B-Alert <sup>+</sup> , BHBT, Stellar, Mobita, BioNomadix/BN-Logger
4.4.1	10.05.2015	10, 8.x & 7	MP150/36R, B-Alert, BHBT, Stellar, Mobita, BioNomadix/BN-Logger
4.4	05.22.2014	8, 7 & Vista	MP150/36R, B-Alert, BHBT, Stellar, Mobita, BioNomadix
4.3.1	08.06.2013	8, 7, &	MP150/36R, B-Alert, BHBT

#### macOS

Version	Build Date	Supported OS	Hardware
5.0.8	09.23.2022	10.13-12*	MP160/150/36R/ BioNomadix/BN-Logger/Smart Center
5.0.6	06.02.2021	10.10-10.15	MP160/150/36R/ BioNomadix/BN-Logger/Smart Center
5.0.5	08.16.2019	10.9-10.14	MP160/150/36R/ BioNomadix/BN-Logger/Smart Center
5.0.4	01.10.2019	10.9-10.14	MP160/150/36R/ BioNomadix/BN-Logger/Smart Center
5.0.3	09.07.2018	10.9-10.13	MP160/150/36R/ BioNomadix/BN-Logger/Smart Center
5.0.2.1	11.27.2017	10.9-10.12	MP160/150/36R/ BioNomadix/BN-Logger/Smart Center
5.0.1	10.31.2016	10.9-10.12	MP160/150/36R/ BioNomadix/BN-Logger
5.0	07.28.2016	10.9-10.11	MP160/150/36R/ BioNomadix/BN-Logger
4.4.2	04.08.2016	10.9-10.11	MP150/36R/ BioNomadix/BN-Logger
4.4.1	10.05.2015	10.7-10.10	MP150/36R/ BioNomadix/BN-Logger
4.4	05.22.2014	10.6-10.9	MP150/36R/ BioNomadix
4.3.1	08.05.2013	10.6-10.8	MP150/36R
4.3	11.05.2012	10.6-10.8	MP150/36R
4.2	05.31.2011	10.5-10.7	MP150/36R
4.1.1	04.22.2010	10.5-10.6	MP150/36R
4.1	07.24.2009	10.5-10.6	MP150/36R
3.9.2	03.14.2007	10.4-10.6	MP150/100/35
3.9.1	08.16.2006	10.3-10.4	MP150/100
3.9.0	11.10.2004	10.2-10.3.9	MP150/100

		Vista	
4.3	11.09.2012	8, 7, & Vista	MP150/36R, B-Alert, BHBT
4.2	05.31.2011	7 & Vista	MP150/36R, B-Alert
4.1.1	04.23.2010	7 & Vista	MP150/36R
4.1BH	04.06.2010	7 & Vista	BioHarness ISM
4.1	07.24.2009	7 & Vista	MP150/36R
4.0	07.22.2008	7 & Vista	MP150
Version	Build Date	Supported OS	Hardware
3.9BH	12.12.2008	Vista & XP	BioHarness ISM
3.9.1	07.11.2007	7, Vista & XP	MP150/100
3.9.0.1	04.20.2007	2000 & XP	MP150/100
3.9.0	11.17.2006	2000 & XP	MP150/100
3.8.2	02.21.2006	2000 & XP	MP150/100
3.8.1**	08.20.2004	2000 & XP	MP150/100

Version	Build Date	Supported OS	Hardware
3.7.0	04.04.2001	8.6-10.1	MP150/100
3.5.0	05.15.2001	8.6-9.1	MP150/100
3.2.7	04.05.2000	7x-9.0	MP100

†B-Alert functionality in AcqKnowledge 4.4.2 is officially supported in Windows 7 & 8 only.

\* Only MP36 units with serial numbers beginning "MP36A" are compatible with macOS 11 and later.

\*\*3.8.1 GLP Version – XP only – **(Discontinued 2006)**

## AcqKnowledge Translations: -CS, -CT, -ES, -FR, -IT, -JP (Windows only) (Discontinued 2008)

Version	Build Date	Supported OS	Hardware
3.8.2-XX	03.29.2007	2000-XP	MP150/MP100



**Biopac  
Student Lab**

English (see page 3 for translations)

### Windows

Version	Build Date	Supported OS	Hardware
4.1.7	11.03.2022	10 & 11	MP36A/36/35/46/41
4.1.6	02.04.2021	10 & 8.x	MP36/35/46/45/41
4.1.5	11.10.2020	10 & 8.x	MP36/35/45/41
4.1.3	02.15.2018	10 & 8.x	MP36/35/45
4.1.2	01.11.2017	10, 8.x & 7	MP36/35/45
4.1.1	10.14.2015	10, 8.x & 7	MP36/35/45
4.1	02.12.2015	8 & 7	MP36/35/45
4.0.3	07.31.2014	8, 7 & Vista	MP36/35/45
4.0.2	01.07.2014	8, 7 & Vista	MP36/35/45
4.0.1	06.26.2013	8, 7, & Vista	MP36/35/45
4.0	04.04.2012	7 & Vista	MP36/35/45
3.7.7	11.25.2010	7, Vista & XP	MP36/35/30/45
3.7.6.1†	03.08.2010	Vista & XP	MP36/35/30
3.7.6†	09.17.2009	Vista & XP	MP36/35/30
3.7.5	06.05.2009	Vista & XP	MP45

### macOS

Version	Build Date	Supported OS**	Hardware
4.1.7	11.04.2022	10.13-12***	MP36A/36/35/46/41
4.1.6	02.04.2021	10.15-10.14	MP36/35/46/45/41
4.1.5	11.10.2020	10.15-10.10	MP36/35/45/41
4.1.3	02.15.2018	10.13-10.10	MP36/35/45
4.1.2	01.12.2017	10.12-10.9	MP36/35/45
4.1.1	10.14.2015	10.9-10.11	MP36/35/45
4.1	02.12.2015	10.7-10.9	MP36/35/45
4.0.3	07.03.2014	10.6-10.9	MP36/35/45
4.0.2	01.07.2014	10.6-10.8	MP36/35/45
4.0.1	06.26.2013	10.6-10.8	MP36/35/45
4.0	04.04.2012	10.5-10.7	MP36/35/45
3.7.6†	03.08.2010	10.4.10.6	MP36/35
3.7.5	06.05.2009	10.4-10.6	MP45

Combo Install			
3.7.5	12.25.2008	Vista & XP	MP45
3.7.3	08.26.2008	Vista & XP	MP36/35/30
3.7.2	05.25.2007	Vista & XP	MP35/30
3.7.1.1	05.16.2006	XP	MP35/30
Version	Build Date	Supported OS	Hardware
3.7.1	04.29.2005	XP, 2000, 98SE	MP35/30
3.7.0	03.24.2004 to 08.10.2004	XP, 2000, 98SE	MP30
3.6.7	08.13.2001	98, ME, 2000	MP30

3.7.3	12.10.2008	10.4-10.6	MP36/35
3.7.2	06.08.2007	10.3-10.4	MP35/30
3.7.1.1	05.15.2006	10.3-10.4	MP35/30
3.7.1	11.01.2005	10.3-10.4	MP35/30
3.6.6	09.16.2002	8.0-9.0	MP30
Version	Build Date	Supported OS**	Hardware
3.0.7	1.20.2002	Web Version	MP30
3.0.7	08.17.2001	7x-8.6	MP30
3.0.6	05.05.2000	7x-8.6 USB	MP30

†Combined BSL 3.7.6 Lessons & BSL PRO 3.7.3

\*\*\* The MP35 is not compatible with macOS versions later than 10.15. Only MP36 units with serial numbers beginning "MP36A" are compatible with macOS 11 and later.



**Biopac  
Student Lab**

**Translations** (see page 2 for English)

Version	Language	Build Date	Supported OS	Hardware
<b>Windows &amp; Mac OS X</b>				
4.1.7	ES Spanish FR French JP Japanese CT Chinese Traditional CS Chinese Simplified	Win 11.03.2022 Mac 11.04.2022	10 & 11 OS X 10.13-12	MP36A/36/35/ 46/41
4.1.6	ES Spanish FR French JP Japanese CT Chinese Traditional CS Chinese Simplified	02.04.2021	10 & 8.x OS X 10.14-10.15	MP36/35/ 46/45/41
4.1.5	ES Spanish FR French JP Japanese CT Chinese Traditional CS Chinese Simplified	11.10.2020	10 & 8.x OS X 10.10-10.15	MP36/35/45/41
4.1.3	ES Spanish FR French JP Japanese CT Chinese Traditional	02.15.2018	10, 8.x, 7 OS X 10.10-10.13	MP36/35/45
4.1.2	ES Spanish FR French JP Japanese	01.11.2017	10, 8 .x & 7 OS X 10.9-10.12	MP36/35/45
4.1.1	ES Spanish FR French JP Japanese	10.14.2015	10, 8.x & 7 OS X 10.9-10.11	MP36/35/45
4.1	ES Spanish	02.12.2015	8 & 7	MP36/35/45

	FR French		OS X 10.7-10.9	
<b>4.0.3</b>	ES Spanish FR French	07.31.2014	8, 7 & Vista OS X 10.6-10.9	MP36/35/45
<b>4.0.2</b>	FR French	01.07.2014	8, 7 & Vista OS X 10.6-10.8	MP36/35/45
<b>Windows only</b>				
<b>3.7.6 Lessons &amp; PRO 3.7.3</b>	RU Russian	08.15.2012	7, Vista, & XP	MP36/35/30
<b>3.7.5</b>	CT Chinese Traditional	06.06.2011	7, Vista & XP	MP45
	CS Chinese Simplified	06.30.2011		
	FR French	12.14.2010		
	ES Spanish	08.30.2011		
	RU Russian	01.06.2012		
<b>Windows only</b>				
<b>3.7.3</b>	JP Japanese	05.07.2010	7, Vista & XP	MP36/35/30
	CT Chinese Traditional	05.26.2010		
	CS Chinese Simplified	06.07.2010		
	IT Italian	06.14.2010		
	FR French	04.19.2010		
	ES Spanish	04.26.2010		
	RU Russian	02.06.2012		
<b>3.7.2</b>	JP Japanese	03.14.2008	Vista & XP	MP35/30
	CT Chinese Traditional	03.13.2008		
	FR French	03.07.2008		
	ES Spanish	05.21.2008		



### Windows

Version	Build Date	Supported OS	Hardware
1.1.1	11.27.2007	10/7/Vista/XP	MP40
1.1	02.08.2006	XP	
1.0	12.22.2005	XP	

### Mac (Discontinued July 2012)

Version	Build Date	Supported OS	Hardware
1.1.1	04.26.2007	10.4-10.5	MP40
1.1	03.09.2006	10.3-10.4	
1.0	12.22.2005	10.3-10.4	